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Today it is more essential than ever for the dentist to take care of the problem of the treatment of the waters.

The main goal to **guarantee** a complete process that is able to **prevent** any disease for patients, providing the highest efficiency and hygiene is becoming a point of success for a dental cabinet.

In addition the dentist is feeling a growing importance in turning the **expensive maintenance** operations into a complete set of automatic ones.

A correct and complete treatment of the waters, when supported by the cooperation of all people working in the cabinet, allows the dentist to:

- **guarantee** the health of patients and of all personnel in the cabinet
- **make** all cleaning and sterilizing operations become **easier**
- **improve** the efficiency of the equipments and autoclaves preventing many functional problems

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What's about biofilm

A **biofilm** is an aggregate of microorganisms in which cells are stuck to each other and/or to a surface.

Biofilms can contain many different types of microorganism, e.g. **bacteria**, archaea, protozoa, fungi and algae.

This sort of structure is a **high efficiency protection** for those kinds of microorganism.

In many cases the **resistance to the antibiotics has been improved up to 1000 times due to the biofilm presence**.

BIOFILM'S DAMAGE

Biofilm, mostly in case of **dental equipments**, can form a very dangerous **bacterial deposit** that can result in **contamination** of the dental water lines and that can **resist** to many kinds of disinfectants.

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The **water** is used constantly on the devices of **every day use in a dental cabinet**.

A. Potable water that supplies the dental unit is often full of **sediments** or even **toxic materials**. Some minerals and limestone can even **seriously damage** many parts of the dental unit (handpieces, syringes, scalers, etc.); but the most important thing is to prevent bacteria can live within the dental unit water lines forming biofilms to protect themselves and that can made them become resistant to disinfectants.

B. Water is also used to clean all the instruments and, if not properly treated, it can also hel the oxidation process due to the calcium carbonate, holding back the tensioattivi of the common cleaning products.

C. Autoclaves must be supplied with distilled or bi-distilled water or with water treated by Osmosis Systems.

In the potable water there can be present cellular deposits and inorganic salts that can feed and protect possible bacterias inside the dental unit water lines.

Some bacterias, like *Helicobacter*, *Pseudomonas*, *Legionella* and *Mycobacterium*, particularly dangerous in how much cause of serious illnesses, can proliferate and grow up, finding a favorable environment inside the dental unit water lines.

The risk of a serious and dangerous cross contamination involves all surgery instruments, the dental unit water lines and can affect all the personnel working in the dental cabinet and their customers.

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What's about limestone

Normally, water contains calcium and magnesium ions that, when the water temperature increases, start to precipitate and form scales and scums.

Limestone can easily damage pipes and other surfaces of the equipments (dental unit water lines, handpieces, scalers, etc..)

A proper treatment of the water can strongly reduce the limestone and sediments in the water that supplies the dental unit and all other machines in the dental cabinet, also strongly reducing the risks of breakdowns and malfunctions and lowering the maintenance costs.

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
Each solution has its own osmotic pressure. Net movement of solvent is from the less-concentrated (hypotonic) to the more-concentrated (hypertonic) solution, which tends to reduce the difference in concentrations and to make them become equal (isotonic).

THE OSMOSIS

Osmosis is the diffusion of water through a semi-permeable membrane. More specifically, it is the movement of water across a semi-permeable membrane from an area of high water potential (low solute concentration) to an area of low water potential (high solute concentration). It is a physical process in which a solvent (the water) moves, without input of energy, across a semi-permeable membrane (permeable to the solvent, but not the solute) separating two solutions of different concentrations.

If we apply a mechanical pressure to the high solute concentration we can obtain a REVERSE OSMOSIS

REVERSE OSMOSIS

A mechanical pressure is applied to an impure solution to force pure water through a semi-permeable membrane. Reverse osmosis is theoretically the most thorough method of large scale water purification available

THE SEMIPERMEABLE MEMBRANE


It is the most important part of the Osmosis process that allows the treatment of water, filtering the water sediments and bacterias to the smallest dimensions and providing pure, clean and pleasant water.

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The reverse osmosis works by using pressure to force a solution through a membrane, retaining the solute on one side and allowing the pure solvent to pass to the other side. This is made when a pressure is applied to the high solute concentrate and is the contrary of the normal osmotic process.

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Dental Unit Water Lines

Problems concerning the treatment of Dental Unit Water Lines has been solved by manufacturers in many different ways:

Independent Reservoirs

This is the easiest solution that is based on independent reservoirs for distilled water that can be treated with some disinfectants (e.g. chlorhexidine).

Anyway it must be underlined that the water provided for the cup and bowl fill unit and for the cuspidor bowl is not treated and in addition this could be a very expensive solution.

Integrated Disinfection Systems

Some manufacturers developed different types of integrated disinfection systems.

Anyway, also in these cases, the water provided for the cup and bowl fill unit and for the cuspidor bowl is not treated.

In these cases, disinfection is made only by contact between the water and the disinfectant and only inside the dental unit.

The Dental Unit Water lines are not provided with treated water.

Reverse Osmosis Filtration

Reverse Osmosis Filtration Systems provide treated pure water to the Dental Unit Water Lines, to the Cup and Bowl Fill unit and to the Cuspidory Bowl, preventing any problems caused by microorganisms, sediments and limestone that can normally be found in non-treated water.

The efficiency of this kind of Systems can be also improved by the addition of Oxygenated Water (Hydrogen Peroxide - H_2O_2) that, if correctly dosed, is a good biocide and is able to destroy possible formations of biofilm inside the Dental Unit Water Lines

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Hygiene by oxygenated water

Different papers showed that bacterial contamination of dental unit waterlines can reach much higher concentrations than allowed for potable water, according to European Community Standards.

Reasons for that are recognized in direct contamination from bacteria present in potable water that supplies dental units or aspirated from the outer surface of dental handpieces.

Data obtained show that a daily disinfection cycle in combination with concentrations of hydrogen peroxide continuously supplied in dental unit waterlines allowed to reduce bacterial contamination under the level of sensitivity of the adopted microbiological assay.

Treatment of water with a Reverse Osmosis System and the addition of Hydrogen Peroxide can be the best and cheapest solution for the dentist to prevent formation of biofilms inside the Dental Unit Water Lines.

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The **Osmodent** machine is a **multi functional** Reverse Osmosis System that can provide treated pure water up to a maximum of **10 dental units** and up to a maximum of **3 autoclaves** at the same time with separated feeding plug from 0 to 15 microsiemens.

The system is provided with an **External reservoir** to guarantee the continuous feeding of pure water without interruptions (capacity **12 litres**).

Osmodent is also equipped with an internal **UV LED Light Bactericide System** and with an internal **Anti-Flooding system**.

Osmodent can be equipped with the **Automatic Hygien System DOSA** (optional). It should be placed in a **technical room**.

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The **Osmo R** machine is a **multi functional** Reverse Osmosis System that can provide treated pure water up to a maximum of **1 dental unit** (water hardness 15 microsiemens) or up to a maximum of **1 autoclave** (water hardness 0 microsiemens).

The system is provided with an **internal reservoir** to guarantee the continuous feeding of pure water without interruptions (capacity **3,5 litres**).

Osmo R is also equipped with an internal **UV Light Bactericide System** and with an internal **Anti-Flooding system** and can be equipped with an **Hygiene System** (optional).

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The **Osmo R3** machine is a **multi functional** Reverse Osmosis System that can provide treated pure water up to a maximum of **3 dental units** (water hardness 15 microsiemens) and up to a maximum of **1 autoclave** (water hardness 0 microsiemens).

The system is provided with an **internal reservoir** to guarantee the continuous feeding of pure water without interruptions (capacity **10 litres**), with an internal **UV Light Bactericide System** and with an internal **Anti-Flooding system**.

Osmo R and can be equipped with an **Hygiene System** (optional).

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The **Osmo A** machine is an extremely compact, versatile and cheap Reverse Osmosis System that is to be connected directly to 1 autoclave, providing a continuous supply of treated pure water (water hardness 0 microsiemens).

The **Osmo A System** can be installed almost anywhere and it can be completed with an external reservoir - capacity 3,5 litres (optional).

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